

Instrument Detergents and Enzymatic Cleaners



NPH Klenz Neutral
Detergent



Enzymatic Cleaners - foam or sponge



Medzyme Enzymatic Cleaner

Instrument Lubricants



Instrument Cleaning Accessories

(Must be disposable, or cleaned and minimally high level disinfected or sterilized between uses.)



Glutaraldehyde High Level Disinfectant & Test Strips



Cidex OPA (orthophthalaldehyde) and test strips



Biological Indicators



Self-contained BI



Process Challenge Device PCD or Test Pack containing BI



Incubator



Log for BI monitoring

Chemical Indicators for Sterilization



Steam




Steam/Ethylene oxide (EO)



Peracetic Acid (Steris System 1)



Bowie-Dick Air Removal Test Packs



Post HLD Rinsing

(Sterile or sub-micron filtered water is **recommended**)




Sterile Water for Irrigation - 500ml
1000ml, 3000ml volumes



Pall-Aquasafe™ Faucet Water
Filter



**Read labels on all cleaning/ disinfection
products. Ensure directions are followed.**



(Sample enzymatic cleaner label)